

(Demo) Product Specification Cat. No. AG050

Description: Rat monoclonal antibody purified from 5 liters tissue culture supernatant.
 Methodology: Protein G
 Clone Designation: YL 1/2
 Buffer: Phosphate Buffered Saline pH 7.4

Lot Number :	11-179
Volume:	50 ml
Protein Concentration:	1.09 mg/ml ¹
Total Yield:	54.5 milligrams
¹ Protein concentration is determined in a modified dye binding assay using Bovine Gamma Globulin (BGG) as a standard	

