

Custom Goat IgG Purification

Description: Custom Goat IgG purification from serum or plasma.
 Animal Designation: 31-080-2, 31-080-4
 Immunogen: Protein G94, HA27
 Buffer: Phosphate Buffered Saline pH 7.4

Lot Number:	33-025-1	33-025-2
Animal:	31-080-2	31-080-4
Immunogen	G94	HA27
Volume of serum or plasma:	1126 ml	1116 ml
Volume of purified product:	1000 ml	1000 ml
Protein Concentration:	33.5 mg/ml ¹	31.0 mg/ml ¹
Total Yield:	33.5 grams	31.0 grams

¹ Protein concentration is determined in a modified dye binding assay using Bovine Gamma Globulin (BGG) as a standard

